
DE Gebrauchsanweisung EN-US Instructions for use FR Mode d'emploi IT Istruzioni d'uso
ES-MX Manual de instrucciones PT Instruções de uso NL Gebruiksaanwijzing DA Brugsan-
visning NO Bruksanvisning SV Bruksanvisning FI Käyttöohje TR Kullanma talimatı EL دليل
الاستخدام KO 사용 설명서 AR تعليمات الاستخدام


WM 68290j 08/2024 DE, EN-US, FR, IT, ES-MX, PT, NL, DA, NO, SV,
FI, TR, EL, KO, AR

CE 0197

 **Manufacturer**
Löwenstein Medical
Technology GmbH + Co.KG
Kronsaalweg 40
22525 Hamburg, Germany
T: +49 40 54702-0
F: +49 40 54702-461
www.loewensteinmedical.com

CARA
Full Face Mask


WM 68290j

LÖWENSTEIN
medical

LÖWENSTEIN
medical

9 Druck-Flow-Kennlinie

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Symbol	Description
	Indicates the product is a medical device
	Follow the instructions for use
	Permitted temperature range for transport and storage
	Use by date
	Keep out of sunlight
	Lot number
	CE symbol (confirms that the product conforms to the applicable European directives/regulations)

11 Warranty

Löwenstein Medical Technology gives the purchaser of a new original Löwenstein Medical Technology product and of a spare part fitted by Löwenstein Medical Technology a limited manufacturer warranty in accordance with the warranty conditions applicable to the product in question and in accordance with the warranty periods from date of purchase listed below. The warranty conditions are available on the manufacturer's website. We will also send you the warranty conditions on request.

Please bear in mind that any claim under warranty and liability shall be void if neither the accessories recommended in the instructions for use nor original spare parts are used.

In the event of a claim under warranty, contact your specialist dealer.

Product	Warranty periods
Masks including accessories	6 months

12 Declaration of Conformity

The manufacturer Löwenstein Medical Technology GmbH + Co. KG (Kronsaalsweg 40, 22525 Hamburg, Germany) hereby declares that the product complies with the relevant provisions of the Medical Device Regulations (EU) 2017/745. The unabridged text of the Declaration of Conformity can be found on the manufacturer's website.

In the EU: As a user and/or patient, you must report any serious incidents occurring in conjunction with the product to the manufacturer and to the responsible authority.

1 Operation

The figures show the following steps for operating the mask:

- 1 Putting on the mask
- 2 Adjusting the mask
- 3 Removing the mask
- 4 Dismantling the mask
- 5 Assembling the mask

For blind and partially-sighted users

An electronic version of the instructions for use is also available on the manufacturer's website.

2 Introduction

2.1 Intended use

The CARA Full Face mask is used for treating sleep apnea and for non-invasive and non-life-sustaining ventilation of patients with respiratory insufficiency. It serves as a connecting element between the patient and the therapy device.

2.2 Contraindications

The mask must not be used for patients weighing < 30 kg.

The mask may not be used in the following situations: Necessity for immediate intubation, loss of consciousness, acute vomiting.

The mask may be used in the following situations only with particular caution: pressure points and acute injuries to the facial skin, skin allergies in the facial area, facial or nasopharyngeal deformities, acute pain in the facial area, cough reflex restricted or absent, claustrophobia, acute nausea.

Mask is not suitable for use with a nebulizer or as an anesthetic mask.

If you are not certain whether one of these situations applies to you, consult your healthcare professional. Please observe the contraindications in the instructions for use of your device.

2.3 Side effects

The following side effects may occur with use of the mask: Nasal congestion, dry nose, dry mouth in the morning, feeling of pressure in the sinuses, irritated conjunctiva, skin rashes, pressure marks on the face, irritating noises when breathing.

If these side effects occur, contact your healthcare professional.

2.4 Clinical benefit

Transfers the therapeutic efficacy of the therapy device to the patient.

3 Safety

⚠ WARNING! Unusually hazardous situation. If you ignore the following instructions, severe, irreversible or fatal injuries may result.

Risk of injury due to damaged mask parts or those under strain!

- ⇒ Perform a visual inspection before every use and after every cleaning operation. Replace mask parts if necessary.
- ⇒ Note useful life (see section entitled "Technical data").

Risk of injury from re-inhaling CO₂!

- ⇒ Only use the mask when therapy is in progress.
- ⇒ Only use the mask within the quoted therapy pressure range.
- ⇒ Patients unable to remove the mask themselves must be monitored by a nurse.
- ⇒ Check before every use that the openings of the anti-asphyxia valve are clear.
- ⇒ Do not close off exhalation systems.

Risk of injury due to the use of oxygen!

Oxygen can accumulate in clothing, bed linen, and hair. Supplying oxygen without a safety device can lead to fire.

- ⇒ Use an oxygen safety valve.
- ⇒ Follow the instructions for use for the oxygen supply system.
- ⇒ Do not smoke.
- ⇒ Avoid naked flames.
- ⇒ Ventilate the room well.
- ⇒ Keep mask free from oil and grease.
- ⇒ Observe the safety instructions in the instructions for use of the device.

Risk of injury due to insufficient supply to the patient!

- ⇒ Activate low pressure/leakage alarms on device.
- ⇒ Use the correct mask size and check that it fits tightly.
- ⇒ Monitor patients with restricted spontaneous breathing.

Risk of injury from inadequate cleaning!

- ⇒ Clean mask parts before using for the first time (see section entitled "Cleaning and reprocessing").
- ⇒ Clean the mask regularly.
- ⇒ In case of a change of patient in the clinical environment: Follow the document *Mask reprocessing instructions* (see section entitled "Reprocessing").

⚠ CAUTION! Hazardous situation. If you ignore these instructions, mild or moderate injuries may result.

Risk of injury due to allergies!

- ⇒ Be aware of the materials used in the mask (see section entitled "Materials"). Only use a mask following consultation with a healthcare professional.
- ⇒ When selecting the detergent, consider potential allergies.

4 Product description

A diagram of the individual parts can be found on the title page.

1	Headgear	5	Rotating sleeve
2	Connecting element	6	Mask body
3	Elbow	7	Headgear clip
4	Anti-asphyxia valve	8	Mask cushion (two mask cushions provided, depending on model)

Compatible devices

In some device combinations, the actual pressure does not correspond to the therapy pressure displayed by the device. Have the device adjusted by a healthcare professional so that the actual pressure in the mask corresponds to the therapy pressure. This adjustment should be made with the type of mask that will be used during therapy.

Exhalation system

The mask has an integrated exhalation system. The exhaled air escapes through a gap.

Anti-asphyxia valve (AAV)

If the device fails, the anti-asphyxia valve opens and the patient breathes ambient air.

Quick-release cord (optional)

The quick-release cord allows the mask to be released quickly and easily in emergency situations (see SOS illustration).

5 Cleaning and reprocessing

5.1 Clean mask

1. Wash your hands before starting cleaning.
2. Dismantle mask (see Figure 4).
3. Clean the mask by hand (max. 30 °C, 1 ml mild detergent to 1 l water) in accordance with the following table:

Mask part	Frequency	Action
All mask parts	Daily	Soak and wash for 15 minutes and clean for 3 minutes using a soft cleaning brush.
Headgear	Weekly	Wash for 15 minutes.

All parts (exception: elbow with anti-asphyxia valve) can be washed in a dishwasher once a week (max. 70 °C, mild dishwashing detergent, max. 90 minutes program duration, top basket, separate rinse).

4. Rinse all parts with clean water.
5. Allow all parts to air-dry.
6. Carry out a visual inspection for cracks and deformations. Replace damaged parts. Discoloration is harmless.
7. Re-assemble mask (see Figure 5).

5.2 Reprocessing (clinical environment)

In the event of a change of patient, follow the instructions in the document entitled *Information on reprocessing*. The document can be found on the manufacturer's website. We will send you this document on request.

5.3 Disposal

Dispose of the mask in domestic waste. In the clinical environment: Dispose of the mask in accordance with hospital regulations.

6 Troubleshooting

Fault	Cause	Action
Pain due to pressure on the face	Mask too tight.	Adjust headgear so that it is looser.
Draft in the eye	Full face mask fits too loosely.	Adjust headgear so that it is tighter.
	Mask does not fit.	Contact your specialist dealer.
Therapy pressure is not reached.	Mask is not correctly adjusted.	Re-adjust mask.
	Mask cushion is damaged.	Replace mask cushion.
	Circuit is damaged.	Check circuit and the correct fit of circuit.
Therapy pressure not reached.	Anti-asphyxia valve defective.	Replace anti-asphyxia valve.

7 Technical data

Dimensions in mm (W x H x D)	
Size S	93 x 145 x 89
Size M	93 x 160 x 90
Size L	93 x 174 x 91
Size XL	96 x 185 x 100
Weight	
Size S	93 g
Size M	97 g
Size L	102 g
Size XL	109 g
Dead space volume	
Size S	180 ml
Size M	219 ml
Size L	244 ml
Size XL	291 ml
Tube connection: Tapered connection to EN ISO 5356-1	Ø 22 mm (male)
Flow resistance at 50 l/min	0.12 hPa
at 100 l/min	0.5 hPa
AAV flow resistance	
Inhalation at 50 l/min	0.6 hPa
Exhalation at 50 l/min	0.8 hPa
Tolerance: ± 0.2 hPa	
AAV switching pressure	
Open	0.5 hPa
Close	2.2 hPa
Service life	5 years
Useful life	Up to 12 months ¹
Therapy pressure	4 hPa - 25 hPa
Quoted two-figure noise emission value according to ISO 4871:	
Sound pressure level	19 dB(A)
Sound power level	27 dB(A)
Uncertainty factor	3 dB(A)
Temperature:	
Operation	+5 °C to +40 °C
Transport and storage	-20 °C to +70 °C
Standards applied	EN ISO 17510: 2020
Product class to MDR (EU) 2017/745	IIa


¹ Useful life depends on cleaning and on the cleaning agents used, on the amount of time worn daily, on therapy pressure, and on individual secretion of sweat. The use of the second mask cushion (optionally included) will not extend the useful life.

8 Materials

No parts of the mask contain latex, PVC (polyvinyl chloride) or DEHP (diethylhexyl phthalate).

Headgear	CO (cotton), PA (polyamide), PU (polyurethane)
Headgear clip	PA (polyamide)
Quick-release cord	PET (polyethylene terephthalate fiber), PA (polyamide)
Quick-release cord clip	PA (polyamide), POM (polyoxymethylene)
Forehead support	PA (polyamide)
Connecting element	PA (polyamide)
Mask body	PA (polyamide)
Mask cushion	SI (silicone)
Elbow	PA (polyamide)
Rotating sleeve	PA (polyamide)
Anti-asphyxia valve	SI (silicone)
Valve safety device	PP (polypropylene)

9 Characteristic pressure/flow curve

10 Markings and symbols

The following markings and symbols may be applied to the device, accessories or packaging.

Symbol	Description
	Manufacturer and, if necessary, date of manufacture
	Unique device identifier (unique product code for medical devices)
	Order number

Douille rotative	PA (polyamide)
Valve anti-asphyxie	SI (silicone)
Dispositif de sécurité de la valve	PP (polypropylène)

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